




Probe CCNE1/CEN19p FISH Probe  
Product No. GC013

**CCNE1**  
Probe Size Approximately 150kb  
Fluorophore TexRed  
Location 19q12  
Also known as CCNE; pCCNE1; cyclin E1

**CEN19p**  
Probe Size Approximately 570kb  
Fluorophore FITC  
Location 19p13.11



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.