




Probe AKT2/CEN19p Dual Color FISH Probe  
Product No. GC014

**AKT2**

Probe Size Approximately 160kb  
Fluorophore TexRed  
Location 19q13.13  
Also known as PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA

**CEN19p**

Probe Size Approximately 570kb  
Fluorophore FITC  
Location 19p12



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.