




Probe AKT2/CEN19p Dual Color FISH Probe
Product No. GC014

AKT2

Probe Size Approximately 160kb
Fluorophore TexRed
Location 19q13.13
Also known as PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA

CEN19p

Probe Size Approximately 570kb
Fluorophore FITC
Location 19p13.11



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.