





Probe MLL/CEN11p Dual Color FISH Probe  
Product No. GC016

**MLL**

Probe Size Approximately 380kb  
Fluorophore TexRed  
Location 11q23.3  
Also known as HRX; MLL; MLL1; TRX1; ALL-1; CXXC7; HTRX1; MLL1A; WDSTS; MLL-AF9; MLL/GAS7; TET1-MLL

**CEN11p**

Probe Size Approximately 630kb  
Fluorophore FITC  
Location 11p11.12



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.