



Probe CHGA-ITPK1/CEN14q Dual Color FISH Probe  
Product No. GC035

CHGA-ITPK1  
Probe Size Approximately 280kb  
Fluorophore TexRed  
Location 14q32  
Also known as CGA- ITRPK1

CEN14q  
Probe Size Approximately 400kb  
Fluorophore FITC  
Location 14q11.2



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.