

Probe AIB1/CEN20q Dual Color FISH Probe  
Product No. GC039

**AIB1**

Probe Size Approximately 290kb

Fluorophore TexRed

Location 20q12



Also known as ACTR; NCOA3; RAC3; SRC3; pCIP; AIB-1; CTG26; SRC-3; CAGH16; KAT13B; TNRC14; TNRC16; TRAM-1; bHLHe42

**CEN20q**

Probe Size Approximately 680kb

Fluorophore FITC

Location 20q11.22



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.