




Probe ERCC1/CEN19p Dual Color FISH Probe  
Product No. GC043

ERCC1  
Probe Size Approximately 240kb  
Fluorophore TexRed  
Location 19q13.32  
Also known as UV20; COFS4; RAD10

CEN19p  
Probe Size Apporoximately 570kb  
Fluorophore FITC  
Location 19p13.11



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.