




Probe TPX2/CEN20p Dual Color FISH Probe
Product No. GC047

TPX2

Probe Size Approximately 150kb
Fluorophore TexRed
Location 20q11.2
Also known as DIL2; p100; DIL-2; HCTP4; FLS353; HCA519; REPP86; C20orf1; G20orf2; GD:C20orf1

CEN20p

Probe Size Approximately 400kb
Fluorophore FITC
Location 20p11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.