

Probe ATP7A/CENXp Dual Color FISH Probe
Product No. GC061


ATP7A

Probe Size Approximately 350kb
Fluorophore TexRed
Location Xq21.1
Also known as MK; MNK; DSMAX; SMAX3



CENXp

Probe Size Approximately 550kb
Fluorophore FITC
Location Xp11.22

ATP7A



CENXp



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.