



Probe LC3/CEN20p Dual Color FISH Probe  
Product No. GC065

LC3

Probe Size Approximately 400kb  
Fluorophore TexRed  
Location 20q11.22  
Also known as MAP1LC3A; LC3A; ATG8E; MAP1ALC3; MAP1BLC3

CEN20p

Probe Size Approximately 400kb  
Fluorophore FITC  
Location 20p11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.