



Probe p21/CEN6p Dual Color FISH Probe
Product No. GC067

p21

Probe Size Approximately 270kb
Fluorophore TexRed
Location 6p21.2
Also known as CDKN1A; CIP1; SDI1; WAF1; CAP20; CDKN1; MDA-6; p21CIP1

CEN6p

Probe Size Approximately 660kb
Fluorophore FITC
Location 6p12.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.