

Probe SQSTM1/CEN5q Dual Color FISH Probe
Product No. GC069

SQSTM1

Probe Size Approximately 240kb
Fluorophore TexRed
Location 5q35.3
Also known as p60; p62; A170; OSIL; PDB3; ZIP3; p62B; FTDALS3

CEN5q

Probe Size Approximately 700kb
Fluorophore FITC
Location 5q11.2



SQSTM1

178.6 178.8 179.0 179.2 179.4 179.6 179.8 (MB)



CEN5q

51.6 51.8 52.0 52.2 52.4 52.6 52.8 (MB)



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.