





Probe SOX2/CEN3q Dual Color FISH Probe  
Product No. GC078

SOX2  
Probe Size Approximately 380kb  
Fluorophore TexRed  
Location 3q26.3-q27  
Also known as ANOP3; MCOPS3

CEN3q  
Probe Size Approximately 500kb  
Fluorophore FITC  
Location 3q12.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.