



Probe RET/CEN10p Dual Color FISH Probe  
Product No. GC082

RET

Probe Size Approximately 200kb  
Fluorophore TexRed  
Location 10q11.2  
Also known as PTC; MTC1; HSCR1; MEN2A; MEN2B; RET51; CDHF12; CDHR16; RET-ELE1

CEN10p

Probe Size Approximately 640kb  
Fluorophore FITC  
Location 10p11.22



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.