

Probe MODY1/CEN20q Dual Color FISH Probe
Product No. GC086

PTPRT

Probe Size Approximately 400kb

Fluorophore TexRed

Location 20q13.12




Also known as TCF; HNF4; MODY; FRTS4; MODY1; NR2A1; TCF14; HNF4a7; HNF4a8; HNF4a9; NR2A21; HNF4alpha

CEN20q

Probe Size Approximately 680kb

Fluorophore FITC

Location 20q11.22



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.