





Probe PDGFRA/CEN4p Dual Color FISH Probe
Product No. GC091

PDGFRA
Probe Size Approximately 240kb
Fluorophore TexRed
Location 4q12
Also known as CD140A; PDGFR2; PDGFR-2; RHEPDGFRA

CEN4p
Probe Size Approximately 560kb
Fluorophore FITC
Location 4p12



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.