

Probe AXL/CEN19q FISH Probe  
Product No. GC093


**AXL**

Probe Size Approximately 290kb  
Fluorophore TexRed  
Location 19q13.2  
Also known as ARK; UFO; JTK11; Tyro7




**CEN19q**

Probe Size Approximately 430kb  
Fluorophore FITC  
Location 19q12

AXL



CEN19q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.