




Probe STK15/CEN20p Dual Color FISH Probe  
Product No. GC113

STK15

Probe Size Approximately 290kb  
Fluorophore TexRed  
Location 20q13  
Also known as AIK; ARK1; AURA; BTAK; STK6; STK7; AURKA; AURORA2; PPP1R47

CEN20p

Probe Size Approximately 400kb  
Fluorophore FITC  
Location 20p11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.