



Probe BRAF1/CEN7q Dual Color FISH Probe  
Product No. GC115

**BRAF1**

Probe Size Approximately 390kb  
Fluorophore TexRed  
Location 7q34  
Also known as NS7; BRAF; RAFB1; B-RAF1

**CEN7q**

Probe Size Approximately 800kb  
Fluorophore FITC  
Location 7q11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.