

Probe MRP1/CEN16q Dual Color FISH Probe  
Product No. GC121


**MRP1**

Probe Size Approximately 410kb  
Fluorophore TexRed  
Location 16q21  
Also known as MRP; ABCC; GS-X; ABCC1; ABC29




**CEN16q**

Probe Size Approximately 700kb  
Fluorophore FITC  
Location 16q12.1

**MRP1**



**CEN16q**



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.