




Probe JAZF1/CEN7p Dual Color FISH Probe(Williams Region)
Product No. GC126

JAZF1

Probe Size Approximately 460kb
Fluorophore TexRed
Location 7p15.2-p15.1
Also known as TIP27; ZNF802

CEN7p

Probe Size Approximately 700kb
Fluorophore FITC
Location 7p12.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.