







Probe MLH1/CEN3q Dual Color FISH Probe  
Product No. GC130

**MLH1**  
Probe Size Approximately 290kb  
Fluorophore TexRed  
Location 3p21.3  
Also known as FCC2; COCA2; HNPCC; hMLH1; HNPCC2

**CEN3q**  
Probe Size Approximately 500kb  
Fluorophore FITC  
Location 3q12.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.