

Probe WT1/CEN11p Dual Color FISH Probe
Product No. GC143




WT1

Probe Size Approximately 170kb
Fluorophore TexRed
Location 11p13
Also known as GUD; AWT1; WAGR; WT33; NPHS4; WIT-2; EWS-WT1

CEN11p

Probe Size Approximately 630kb
Fluorophore FITC
Location 11p11.12

WT1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.