




Probe PTP4A3/CEN8q Dual Color FISH Probe  
Product No. GC144

**PTP4A3**

Probe Size Approximately 140kb  
Fluorophore TexRed  
Location 8q24.3  
Also known as PRL3; PRL-3; PRL-R

**CEN8q**

Probe Size Approximately 520kb  
Fluorophore FITC  
Location 8q11.21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.