

Probe SRC/CEN20p Dual Color FISH Probe  
Product No. GC149





SRC

Probe Size Approximately 150kb  
Fluorophore TexRed  
Location 20q12-q13  
Also known as ASV; SRC1; c-SRC; p60-Src

CEN20p

Probe Size Approximately 400kb  
Fluorophore FITC  
Location 20p11.21

SRC



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.