





Probe NTRK2/CEN9q Dual Color FISH Probe  
Product No. GC154

**NTRK2**  
Probe Size Approximately 460kb  
Fluorophore TexRed  
Location 9q22.1  
Also known as TRKB; trk-B; GP145-TrkB

**CEN9q**  
Probe Size Approximately 470kb  
Fluorophore FITC  
Location 9q21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.