



Probe AF4/CEN4p Dual Color FISH Probe  
Product No. GC163

AF4  
Probe Size Approximately 450kb  
Fluorophore TexRed  
Location 4q21  
Also known as AFF1; PBM1; MLLT2



CEN4p  
Probe Size Approximately 560kb  
Fluorophore FITC  
Location 4p13



~450kb



~560kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.