

Probe FLT3/CEN13q Dual Color FISH Probe
Product No. GC170

FLT3

Probe Size Approximately 260kb
Fluorophore TexRed
Location 13q12
Also known as FLK2; STK1; CD135; FLK-2



CEN13q

Probe Size Approximately 550kb
Fluorophore FITC
Location 13q12.11

FLT3



CEN13q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.