


Probe DDR2/CEN1p Dual Color FISH Probe
Product No. GC176

DDR2
Probe Size Approximately 260kb
Fluorophore TexRed
Location 1q23.3
Also known as TKT; MIG20a; NTRKR3; TYRO10



CEN1p
Probe Size Approximately 780kb
Fluorophore FITC
Location 1p13.3



~260kb



~780kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.