

Probe MIIP/CEN1p Dual Color FISH Probe  
Product No. GC177


**MIIP**

Probe Size Approximately 170kb  
Fluorophore TexRed  
Location 1p36.22  
Also known as IIP45





**CEN1p**

Probe Size Approximately 780kb  
Fluorophore FITC  
Location 1p13.3

MIIP



CEN1p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.