

Probe ATRX/CENXp Dual Color FISH Probe  
Product No. GC202

ATRX

Probe Size Approximately 300kb

Fluorophore TexRed

Location Xq21.1

Also known as JMS; SHS; XH2; XNP; ATR2; SFM1; MRX52; RAD54; MRXHF1; RAD54L; ZNF-HX


CENXp

Probe Size Approximately 550kb




Fluorophore FITC

Location Xp11.22

ATRX



CENXp



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.