

Probe MTAP/CEN9q Dual Color FISH Probe  
Product No. GC206


**MTAP**

Probe Size Approximately 210kb  
Fluorophore TexRed  
Location 9p21.3  
Also known as BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249




**CEN9q**

Probe Size Approximately 470kb  
Fluorophore FITC  
Location 9q21

MTAP



CEN9q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.