




Probe RREB1/CEN6p Dual Color FISH Probe
Product No. GC230

RREB1

Probe Size Approximately 260kb
Fluorophore TexRed
Location 6p24.3
Also known as HNT; FINB; LZ321; Zep-1; RREB-1

CEN6p

Probe Size Approximately 660kb
Fluorophore FITC
Location 6p12.1



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.