




Probe SOD2 /CEN6p Dual Color FISH Probe  
Product No. GC108

SOD2  
Probe Size Approximately 180kb  
Fluorophore TexRed  
Location 6q25.3  
Also known as IPOB; MNSOD; MVCD6

CEN6p  
Probe Size Approximately 660kb  
Fluorophore FITC  
Location 6p12.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.