



Probe CD273/CEN9q Dual FISH Probe
Product No. GC198

CD273

Probe Size Approximately 130kb
Fluorophore TexRed
Location 9p24
Also known as B7DC; Btdc; PDL2; PDCD1LG2; PD-L2; PDCD1L2; bA574F11.2

CEN9q

Probe Size Approximately 470kb
Fluorophore FITC
Location 9q21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.