




Probe c-met/CEN7 Dual Color FISH Probe  
Product No. GC215

c-met  
Probe Size Approximately 270kb  
Fluorophore TexRed  
Location 7q31  
Also known as HGFR; AUTS9; RCCP2; MET; DFNB97

CEN7(D7Z1)  
Probe Size Approximately 3Mb  
Fluorophore FITC  
Location 7p11.1-q11.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.