




Probe EGFR/CEN7 Dual Color FISH Probe
Product No. GC216

EGFR
Probe Size Approximately 320kb
Fluorophore TexRed
Location 7p12
Also known as ERBB; HER1; mENA; ERBB1; PIG61; NISBD2

CEN7(D7Z1)
Probe Size Approximately 3Mb
Fluorophore FITC
Location 7p11.1-q11.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.