

Probe 3-Color FISH Probe (CEN13q, CEN18q, CEN21q)  
Product No. CM3002

CEN13q

Probe Size Approximately 550kb  
Fluorophore FITC  
Location 13q12.11


CEN18q

Probe Size Approximately 660kb  
Fluorophore TexRed  
Location 18q11.2


CEN21q (DYRK1A Gene /Down syndrome critical region)

Probe Size Approximately 550kb  
Fluorophore R6G  
Location 21q22.13

CEN13q



CEN18q



CEN21q (DYRK1A Gene /Down syndrome critical region)

