

Probe 5-Color FISH Probe (CEN8q, CEN14q, CEN15q, CENXp, CENYq)  
Product No. CM5003

CEN8q

Probe Size Approximately 520kb  
Fluorophore TexRed  
Location 8q11.21

CEN14q

Probe Size Approximately 400kb  
Fluorophore DEAC  
Location 14q11.2

CEN15q

Probe Size Approximately 680kb  
Fluorophore Cy5  
Location 15q11.2

CENXp

Probe Size Approximately 550kb  
Fluorophore FITC  
Location Xp11.23

CENYq

Probe Size Approximately 470kb  
Fluorophore R6G  
Location Yq11.221



