

Probe 5-Color FISH Probe (CEN8q, CEN14q, CEN15q, CENXp, CENYq)
Product No. CM5003

CEN8q

Probe Size Approximately 520kb
Fluorophore TexRed
Location 8q11.21

CEN14q

Probe Size Approximately 400kb
Fluorophore DEAC
Location 14q11.2

CEN15q

Probe Size Approximately 680kb
Fluorophore Cy5
Location 15q11.2

CENXp

Probe Size Approximately 550kb
Fluorophore FITC
Location Xp11.23

CENYq

Probe Size Approximately 470kb
Fluorophore R6G
Location Yq11.221



