

Probe 3-Color FISH Probe (CEN21q, CENXp, CENYq)
Product No. CM3003

CEN21q (DYRK1A Gene /Down syndrome critical region)

Probe Size Approximately 550kb
Fluorophore R6G
Location 21q22.13


CENXp

Probe Size Approximately 550kb
Fluorophore FITC
Location Xp11.23


CENYq

Probe Size Approximately 470kb
Fluorophore TexRed
Location Yq11.221

CEN21q (DYRK1A Gene /Down syndrome critical region)



CENXp



CENYq

