



Probe BRAF1 Split Dual Color FISH Probe
Product No. SP060

BRAF1
Fluorophore TexRed, FITC
Location 7q34
Also known as NS7; B-raf; BRAF1; RAFB1; B-RAF1

Storage 4°C



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.