

Probe IRF4, DUSP22 Split Dual Color FISH Probe  
Product No. SP079



IRF4, DUSP22

Fluorophore TexRed, FITC

Location 6p25.3

Also known as MUM1; LSIRF; SHEP8; NF-EM5

VHX; JKAP; JSP1; MKPX; JSP-1; MKP-x; LMWDSP2; LMW-DSP2



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.