





Probe BRAF1/LMO7, DY Translocation Dual Color FISH Probe
Product No. TR028

BRAF1
Probe Size Approximately 1050kb
Fluorophore FITC
Location 7q34
Also known as NS7; B-raf; BRAF; RAFB1; B-RAF1

LMO7
Probe Size Approximately 870kb
Fluorophore TexRed
Location 13q22.2
Also known as LOMP; FBX20; FBXO20

Storage 4°C



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.