





Probe RREB1/MKL2, SY Translocation FISH Probe  
Product No. TR043

**RREB1**  
Probe Size Approximately 610kb  
Fluorophore FITC  
Location 6p24.3  
Also known as HNT; FINB; LZ321; Zep-1; RREB-1

**MKL2**  
Probe Size Approximately 560kb  
Fluorophore TexRed  
Location 16p13.12  
Also known as MRTFB; MRTF-B; NPD001



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.