




Probe TFEB/MALAT1, SY Dual Color Translocation FISH Probe
Product No. TR046

TFEB
Probe Size Approximately 470 kb
Fluorophore FITC
Location 6p21.2
Also known as TCFEB; BHLHE35; ALPHATFEB

MALAT1
Probe Size Approximately 640 kb
Fluorophore TexRed
Location 11q13.1
Also known as HCN; NEAT2; PRO2853; LINC00047; NCRNA00047



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.