

Probe NPM1/RARA,DY Translocation Dual Color FISH Probe  
Product No. TR018


RARA

Probe Size Approximately 1200kb  
Fluorophore FITC  
Location 17q21.2  
Also known as RAR; NR1B1




NPM1

Probe Size Approximately 1010kb  
Fluorophore TexRed  
Location 5q35.1  
Also known as B23; NPM

RARA



NPM1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.