

Probe VEGFA/CEN6p Dual Color FISH Probe
Product No. GC128


VEGFA

Probe Size Approximately 240kb
Fluorophore TexRed
Location 6p21.1
Also known as VPF; VEGF; MVCD1




CEN6p

Probe Size Approximately 660kb
Fluorophore FITC
Location 6p12.1

VEGFA



CEN6p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.