

Probe GAPDH/CEN12q FISH Probe  
Product No. GC253

GAPDH  
Probe Size Approximately 530kb  
Fluorophore TexRed  
Location 12p13.31  
Also known as G3PD; GAPD; HEL-S-162eP

CEN12q  
Probe Size Approximately 580kb  
Fluorophore FITC  
Location 12q12

GAPDH  
102.0 102.2 102.4 102.6 102.8 103.0 103.2 (MB)





~530kb

CEN12q  
43.2 43.4 43.6 43.8 44.0 44.2 44.4 (MB)



~580kb



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.