

Probe LAMC2/CEN1q Dual Color FISH Probe
Product No. GC259


LAMC2

Probe Size Approximately 300kb
Fluorophore TexRed
Location 1q25.3
Also known as B2T; CSF; EBR2; BM600; EBR2A; JEB3A; JEB3B; LAMB2T; LAMNB2




CEN1q

Probe Size Approximately 600kb
Fluorophore FITC
Location 1q21.3

LAMC2



CEN1q



representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus